

DATA REPORT - SVP1 AND AMPOULE LINE EMS

Date Time	Temperature	RH	Diffrential Pressure	Diffrential Pressure	Diffrential Pressure
DD-MM-YYYY HH:MM:SS	Washing and Sterilization Area (NMT 30°C)	Washing and Sterilization Area (NMT 80%)	Washing and Sterilization Area w.r.t Unit Preparation Room (NLT 6 Pa)	Washing and Sterilization Area w.r.t Ampoule Decartoning Room (NLT 12 Pa)	Washing and Sterilization Area w.r.t Injectable main corridor (NLT 12 Pa)
25-05-2022 16:58:00	17.6	76	20	16	37
25-05-2022 17:08:00	17.6	76.1	15	15	37
25-05-2022 17:18:00	17.8	80.2	20	16	37
25-05-2022 17:28:00	18.1	80.5	21	15	38
25-05-2022 17:38:00	18.2	80.7	21	16	37
25-05-2022 17:48:00	18.2	81	20	16	37
25-05-2022 17:58:00	18	81	20	14	33
25-05-2022 18:08:00	18	80.5	20	16	37
25-05-2022 18:18:00	18	77	21	14	33
25-05-2022 18:28:00	17.9	76.5	21	16	37
25-05-2022 18:38:00	17.9	76.4	22	16	37
25-05-2022 18:48:00	17.8	76.4	22	16	37
25-05-2022 18:58:00	17.8	76.4	22	16	37
25-05-2022 19:08:00	17.8	76.4	22	16	37
25-05-2022 19:18:00	17.8	76.6	22	16	36
25-05-2022 19:28:00	17.9	76.6	21	16	37
25-05-2022 19:38:00	17.8	76.4	21	16	37
25-05-2022 19:48:00	17.8	77	16	16	37
25-05-2022 19:58:00	18	81	20	14	34
25-05-2022 20:08:00	18.4	79.5	20	15	35
25-05-2022 20:18:00	18.3	77.8	20	16	37
25-05-2022 20:28:00	18.2	77.6	21	16	37
25-05-2022 20:38:00	18.1	77.3	19	15	38
25-05-2022 20:48:00	18.1	77.4	19	16	37

ACTIVITY	NAME	SIGNATURE AND DATE
REVIEWED BY USER		
REVIEWED BY QA		



DATA REPORT - SVP1 AND AMPOULE LINE EMS

Date Time	Temperature	RH	Diffrential Pressure	Diffrential Pressure	Diffrential Pressure
DD-MM-YYYY HH:MM:SS	Washing and Sterilization Area (NMT 30°C)	Washing and Sterilization Area (NMT 80%)	Washing and Sterilization Area w.r.t Unit Preparation Room (NLT 6 Pa)	Washing and Sterilization Area w.r.t Ampoule Decartoning Room (NLT 12 Pa)	Washing and Sterilization Area w.r.t Injectable main corridor (NLT 12 Pa)
25-05-2022 21:28:00	18.3	77.3	0	15	39
25-05-2022 21:38:00	18.3	76.7	0	15	37
25-05-2022 21:48:00	18.2	76.4	19	16	39
25-05-2022 21:58:00	18.3	76.3	19	16	39
25-05-2022 22:08:00	18.4	76.1	19	15	36
25-05-2022 22:18:00	18.4	76	18	16	39
25-05-2022 22:28:00	18.5	75.7	18	16	39
25-05-2022 22:38:00	18.4	75.8	19	16	39
25-05-2022 22:48:00	18.3	77.2	24	17	40
25-05-2022 22:58:00	18.9	76	21	15	36
25-05-2022 23:08:00	18.6	76.5	21	16	39
25-05-2022 23:18:00	18.3	76.8	21	17	41
25-05-2022 23:28:00	18.2	76.8	19	16	40
25-05-2022 23:38:00	18.2	77.8	19	17	42
25-05-2022 23:48:00	18.3	77.6	19	17	41
25-05-2022 23:58:00	18.3	77.1	19	17	41
26-05-2022 00:08:00	18.3	77.4	20	17	41
26-05-2022 00:18:00	18.3	77.2	20	18	41
26-05-2022 00:28:00	18.3	79.2	18	17	41
26-05-2022 00:38:00	18.9	81.2	21	17	41
26-05-2022 00:48:00	19	82.2	22	18	42
26-05-2022 00:58:00	18.9	81	20	16	40
26-05-2022 01:08:00	18.9	81.5	21	15	34
26-05-2022 01:18:00	18.9	82.7	15	17	41

ACTIVITY	NAME	SIGNATURE AND DATE
REVIEWED BY USER		
REVIEWED BY QA		



DATA REPORT - SVP1 AND AMPOULE LINE EMS

Date Time	Temperature	RH	Diffrential Pressure	Diffrential Pressure	Diffrential Pressure
DD-MM-YYYY HH:MM:SS	Washing and Sterilization Area (NMT 30°C)	Washing and Sterilization Area (NMT 80%)	Washing and Sterilization Area w.r.t Unit Preparation Room (NLT 6 Pa)	Washing and Sterilization Area w.r.t Ampoule Decartoning Room (NLT 12 Pa)	Washing and Sterilization Area w.r.t Injectable main corridor (NLT 12 Pa)
26-05-2022 01:58:00	18.7	78.2	19	17	43
26-05-2022 02:08:00	18.7	76.4	20	17	40
26-05-2022 02:18:00	18.5	76.3	19	17	40
26-05-2022 02:28:00	18.4	77.5	21	17	41
26-05-2022 02:38:00	18.4	76.4	19	17	41
26-05-2022 02:48:00	18.4	77.1	18	17	41
26-05-2022 02:58:00	18.4	76.6	19	17	41
26-05-2022 03:08:00	18.4	76.9	19	17	41
26-05-2022 03:18:00	18.4	77.2	19	16	40
26-05-2022 03:28:00	18.4	76.8	20	17	40
26-05-2022 03:38:00	18.4	77.2	19	17	40
26-05-2022 03:48:00	18.5	78.3	0	15	38
26-05-2022 03:58:00	18.6	77.8	20	17	40
26-05-2022 04:08:00	18.7	78.6	20	17	41
26-05-2022 04:18:00	18.8	79.3	19	17	40
26-05-2022 04:28:00	18.9	78.6	20	17	40
26-05-2022 04:38:00	18.9	78.1	20	17	40
26-05-2022 04:48:00	19	76	20	17	40
26-05-2022 04:58:00	18.5	76.9	20	17	40
26-05-2022 05:08:00	18.4	77.3	21	17	41
26-05-2022 05:18:00	18.3	77.4	21	17	40
26-05-2022 05:28:00	18.2	77.5	21	17	40
26-05-2022 05:38:00	18.3	77.4	21	17	41
26-05-2022 05:48:00	18.2	77	21	17	40

ACTIVITY	NAME	SIGNATURE AND DATE
REVIEWED BY USER		
REVIEWED BY QA		



DATA REPORT - SVP1 AND AMPOULE LINE EMS

Date Time	Temperature	RH	Diffrential Pressure	Diffrential Pressure	Diffrential Pressure
DD-MM-YYYY HH:MM:SS	Washing and Sterilization Area (NMT 30°C)	Washing and Sterilization Area (NMT 80%)	Washing and Sterilization Area w.r.t Unit Preparation Room (NLT 6 Pa)	Washing and Sterilization Area w.r.t Ampoule Decartoning Room (NLT 12 Pa)	Washing and Sterilization Area w.r.t Injectable main corridor (NLT 12 Pa)
26-05-2022 06:28:00	18.5	79	19	10	32
26-05-2022 06:38:00	18.7	79.6	19	16	37
26-05-2022 06:48:00	18.8	79.1	20	17	38
26-05-2022 06:58:00	18.9	84.3	19	17	41
26-05-2022 07:08:00	18.9	80.8	19	17	40
26-05-2022 07:18:00	18.8	80.3	18	17	40
26-05-2022 07:28:00	18.8	79.1	20	15	35
26-05-2022 07:38:00	19.2	79.2	20	19	45
26-05-2022 07:48:00	19.4	75.3	18	16	39
26-05-2022 07:58:00	19.1	75.3	18	16	39
26-05-2022 08:08:00	18.9	75.6	19	16	39
26-05-2022 08:18:00	18.8	75.9	19	15	35
26-05-2022 08:28:00	18.7	76.1	18	17	39
26-05-2022 08:38:00	18.6	76.2	18	16	38
26-05-2022 08:48:00	18.6	76.7	19	15	36
26-05-2022 08:58:00	18.5	76.9	19	15	35
26-05-2022 09:08:00	18.5	76.4	19	16	39
26-05-2022 09:18:00	18.4	76.8	19	16	38
26-05-2022 09:28:00	18.4	77.1	19	16	39
26-05-2022 09:38:00	18.4	78.6	19	14	35
26-05-2022 09:48:00	18.7	79.4	19	16	39
26-05-2022 09:58:00	18.7	82	19	17	39
26-05-2022 10:08:00	18.8	82.3	18	16	38
26-05-2022 10:18:00	19	88.5	18	16	39

ACTIVITY	NAME	SIGNATURE AND DATE
REVIEWED BY USER		
REVIEWED BY QA		



DATA REPORT - SVP1 AND AMPOULE LINE EMS

Date Time	Temperature	RH	Diffrential Pressure	Diffrential Pressure	Diffrential Pressure
DD-MM-YYYY HH:MM:SS	Washing and Sterilization Area (NMT 30°C)	Washing and Sterilization Area (NMT 80%)	Washing and Sterilization Area w.r.t Unit Preparation Room (NLT 6 Pa)	Washing and Sterilization Area w.r.t Ampoule Decartoning Room (NLT 12 Pa)	Washing and Sterilization Area w.r.t Injectable main corridor (NLT 12 Pa)
26-05-2022 10:58:00	20.9	84.7	20	17	40
26-05-2022 11:08:00	21.1	82.9	21	17	40
26-05-2022 11:18:00	21.3	81.1	20	17	40
26-05-2022 11:28:00	21.4	80.1	16	17	41
26-05-2022 11:38:00	21.6	78.5	19	17	40
26-05-2022 11:48:00	21.7	77.2	20	16	39
26-05-2022 11:58:00	21.8	76.3	20	17	40
26-05-2022 12:08:00	22.4	74.7	19	17	41
26-05-2022 12:18:00	22.7	72.5	20	16	40
26-05-2022 12:28:00	22.6	70.7	20	15	36
26-05-2022 12:38:00	22.3	71.1	20	17	40
26-05-2022 12:48:00	22.4	71.4	20	17	41
26-05-2022 12:58:00	22.4	70.6	20	17	40
26-05-2022 13:08:00	22.5	70.7	20	15	37
26-05-2022 13:18:00	22.5	70	20	17	40
26-05-2022 13:28:00	22.6	70.3	20	17	41
26-05-2022 13:38:00	22.7	72.9	20	17	41
26-05-2022 13:48:00	22.7	70.6	19	16	37
26-05-2022 13:58:00	22.7	70.4	19	16	37
26-05-2022 14:08:00	22.8	70.1	20	17	41
26-05-2022 14:18:00	22.8	69.7	19	15	36
26-05-2022 14:28:00	22.9	68.4	36	18	41
26-05-2022 14:38:00	23	74.6	20	17	39
26-05-2022 14:48:00	23.1	78	0	14	36

ACTIVITY	NAME	SIGNATURE AND DATE
REVIEWED BY USER		
REVIEWED BY QA		



DATA REPORT - SVP1 AND AMPOULE LINE EMS

WASHING AND STERILIZATION AREA

Date Time	Temperature	RH	Diffrential Pressure	Diffrential Pressure	Diffrential Pressure
DD-MM-YYYY HH:MM:SS	Washing and Sterilization Area (NMT 30°C)	Washing and Sterilization Area (NMT 80%)	Washing and Sterilization Area w.r.t Unit Preparation Room (NLT 6 Pa)	Washing and Sterilization Area w.r.t Ampoule Decartoning Room (NLT 12 Pa)	Washing and Sterilization Area w.r.t Injectable main corridor (NLT 12 Pa)
26-05-2022 15:28:00	23.1	77.4	22	17	41
26-05-2022 15:38:00	23.1	77.3	21	17	41
26-05-2022 15:48:00	23.1	77.7	21	17	41
26-05-2022 15:58:00	23.2	75.9	20	14	36
26-05-2022 16:08:00	23.1	74.2	19	17	40
26-05-2022 16:18:00	23.2	77.3	19	16	41
26-05-2022 16:28:00	23.1	74.1	20	15	38
Min Value	17.6	68.4	0	10	32
Max Value	23.2	88.5	36	19	45

End Of Report: From 25-05-2022 16:28:33 To 26-05-2022 16:28:33

ACTIVITY	NAME	SIGNATURE AND DATE
REVIEWED BY USER		
REVIEWED BY QA		